

BIOSAFETY MANUAL



		Document no.: RTPC2/SOP10	Edition/Revision: 1/0
			1 of 2
Waste Management			Effective Date: March 2013 Date of Revision: NA
Prepared by: Mohamad Supian Mohd Noor (SOP Committee)		Reviewed and approved by: Dr Goh Hoe Han (Head of SOP Committee)	
<u>Purpose of procedure</u> This procedure has been developed to ensure the proper disposal of wastes			
<u>Scope</u> This procedure only applied to the RTPC2 facility management			
<u>Responsibilities</u> Facility technician to ensure proper management of waste of RTPC2 facility			
<u>Steps</u> <ul style="list-style-type: none">• Collect and autoclave waste at least once daily when required• Autoclave or chemically disinfect contaminated protective clothing and containers before reuse for 30 min at 121°C (Operation of autoclave refer to SOP 9)• Autoclaved waste to be transferred out to designated bin for monthly collection by authorised waste collector			
Plant Materials: <ol style="list-style-type: none">1. Cut big plants into smaller manageable pieces2. Place into biohazard autoclave disposable bag and sealed3. Place sealed waste bag into temporary waste disposal bin4. Mark waste bag with autoclave tape before autoclaving5. Steam-sterilize			
Soil (to reuse): <ol style="list-style-type: none">1. Collect soil into biohazard autoclave disposable bag, sealed and mark with autoclave tape2. Steam-sterilized3. Store at an open site to dry			
Culture Materials, Microbial Materials, Liquid Wastes: <ol style="list-style-type: none">1. Cell culture wastes must be decontaminated2. Add household bleach to the container of liquid waste so that the final volume of the solution is 10% bleach.3. Mix gently and allow a minimum of 30 minutes contact time before discarding into the sanitary sewer.			
Pesticide Waste: <ol style="list-style-type: none">1. Label the waste container appropriately according to their chemical class2. Mark the container with the date the wastes were created3. Disposal of pesticide waste shall be conducted by contractor authorized by UKM			

Note:**Verification of effectiveness (monthly)****Seed germination test (tested only for the solid waste programme) :**

Autoclaved soil will be sampled for seed germination test in a pot every month

Print cycle process:

Print detailed cycle process monthly for record

Training required

1. Laboratory induction
2. Emergency response plans
3. Correct use of personal protective equipment and tools

Location

Greenhouse, preparation room and autoclave room