Guideline for specimen collection and transportation for *Helicobacter pylori* culture and susceptibility testing (UKM Medical Centre)

Sample Collection	
p	
Type of specimen	Gastric biopsy (antrum and/or corpus site).
Optimal time of collection	 Before antimicrobial therapy where possible. For patients who are taken antibiotics or <i>H. pylori</i> eradication regimen, collection of specimen is done at least 6 weeks after the therapy has been stopped
Quantity	Depends at the discretion of endoscopist (to increase the yield/sensitivity of culture, several biopsies may required)
Specimen transport & Storage	
Transport medium	Use <i>H. pylori</i> transport medium - The transport medium is available and given upon request - Unused medium should be kept at 4°C (stability maintain for 6 months)
Transport the specimen to: Microbiology Laboratory, 10 th Floor, Preclinical Block, Faculty of Medicine, UKM Medical Centre Jalan Yaacob Latif, Bandar Tun Razak, 56000 Cheras, Kuala Lumpur. Tel: 03-9145 8656 / 8654 / 5920 Fax: 03-9145 9524	Specimen in the transport medium should always be kept at 4°C (in ice) during transportation Send to the lab as soon as possible after collection (preferably within 6 hours) If transportation delays (more than 6 hours): Specimen in the transport medium should be kept at -70°C or -80°C (if available) or -30°C
Operation Hours	Monday to Friday 8.00 am – 4.00pm (Closed : Saturday, Sunday and Public Holidays)
Turn Around Time	30 days
Charge Rate (per sample)	UKM/PPUKM: RM 200.00 Government: RM 200.00 Private: RM 300.00
Methods of reporting result	Hardcopy/ Facsimile / Email / Mail